
# Cole-Parmer®

## GT-200 Series Gelation Timer

Instruction Manual  
6106029(ES)-CPB Version 2.2

**Cole-Parmer®**  
essentials

## CONTENTS


|  | PAGE      |
|--|-----------|
| <b>INTRODUCTION</b>                                | <b>3</b>  |
| <b>GENERAL DESCRIPTION</b>                         | <b>3</b>  |
| <b>PRINCIPLES OF OPERATION</b>                     | <b>4</b>  |
| <b>SAFETY AND INSTALLATION</b>                     | <b>5</b>  |
| Electrical installation                            | 5         |
| <b>SÉCURITÉ ET CONSIGNES D'INSTALLATION</b>        | <b>6</b>  |
| <b>SICHERHEITS- UND INSTALLATIONSINFORMATIONEN</b> | <b>9</b>  |
| <b>INFORMACIÓN DE SEGURIDAD E INSTALACIÓN</b>      | <b>11</b> |
| <b>SPECIFICATION</b>                               | <b>13</b> |
| <b>OPERATION</b>                                   | <b>14</b> |
| Sample preparation                                 | 14        |
| Starting gelation                                  | 14        |
| Accidental tripping during gelation measurement    | 15        |
| Stopping during gelation measurement               | 15        |
| <b>MAINTENANCE, SERVICING &amp; REPAIR</b>         | <b>16</b> |
| Repairs  | 16        |
| Calibration  | 16        |
| <b>ACCESSORIES</b>                                 | <b>17</b> |
| Disposable plungers                                | 17        |
| Stainless steel plungers                           | 17        |
| <b>DECLARATION OF CONFORMITY</b>                   | <b>19</b> |

## INTRODUCTION

**For you own safety and that of others please read and understand all the information in this user manual before using the Gelation Timer.**

## GENERAL DESCRIPTION

The Cole-Parmer Gelation Timer is the ideal instrument for accurate measurement and quality control of gelation of resin or adhesive based samples in the laboratory. All models have digital timers, with the output time in 1 minute increments up to 9999 minutes, for the 1 RPM model (GT-200-5), and in tenths of a minute up to 999.9 minutes, for the 10 RPM model (GT-200-6). The accuracy of the gelation time is  $\pm 2\%$  at a specific temperature. Before leaving the factory each instrument is calibrated and certificated.

## PRINCIPLES OF OPERATION

As polymerisation proceeds in a liquid polymer a point will be reached at which on average there will be one cross link per molecule and at this point gelation will be observed. The gel formed may be soft but it will have the property of a solid of elastic recovery from deformation. This change from liquid to a solid is preceded by a considerable rise in viscosity.

An instrument to measure the time of gelation must be able to discriminate between a state of high viscosity and that of elasticity. A gelation timer operated primarily by change in viscosity would be wrong in principle and, although it might give useful control tests on a particular product, it could not be expected to give times in accordance with the fundamental theoretical analysis of gelation.

The design of the Cole-Parmer Gelation Timer<sup>1</sup> was based on the above considerations. A flat weighted disc or paddle is connected by a link with end play in it to a crank turned by a synchronous motor that falls under gravity in the polymer liquid. The disc is pulled up on the upstroke by the synchronous motor. At the gel point the rigidity of the polymer is sufficient to support the weight of the disc, causing the link to be compressed and to close an electric circuit. This operates a relay which stops the synchronous motor and lights a neon lamp. The gelation time can be read on a digital counter connected to the synchronous motor.

The motion of the disc is resisted both by elastic and by viscous forces but the size, weight and time have been chosen so that the magnitude of a purely viscous force sufficient to trip the mechanism would be far greater than is observed in practice before gelation has taken place.

The work of B.A. Hills<sup>2</sup> has shown that gelation times measured by this instrument with suitable choice of disc diameter, are in excellent agreement with predicted values based on the analysis of P. J. Flory<sup>3</sup>.

Measurements of gelation time with the Cole-Parmer Gelation Timer are thus a very simple method of making fundamental control tests on a polymeric products.

The Cole-Parmer Gelation Timer conforms to the requirements of BS 3532 : 1990 and BS 2782 : Part 8 : Method 835C : 1980 for gelation timing of unsaturated polyester resin systems. It is equally useful for epoxies, urea formaldehyde, resorcinol formaldehyde and other cross linking products.

- 1 N.A. de Bruyne, modern plastics 27, No.9, May, 1950.
- 2 B.A. Hills, J. Oil & Colour Chemists Assoc. 45, 251-260 (1962).
- 3 P.J. Flory, J. American Chemical Soc. (1941) 63, 3083.

The Gelation Timer is designed so that tests may be carried out at constant temperature. The sample cup may stand directly in a laboratory water bath and the gelation timer is relatively narrow so that several instruments may, if necessary, stand together over quite a small bath.

An alternative way of keeping sample cups at constant temperatures is to use a Cole-Parmer Block Heater which is held at accurately controlled temperatures and which has a hole drilled in it of a diameter to accept the sample cup.

Cole-Parmer manufactures a wide range of both water baths and Block heaters suitable for gelation tests at constant temperatures. Details of these are available on request.

## SAFETY AND INSTALLATION

Before discarding the packaging check that all parts are present and correct.

If the equipment is not used in the manner described in this manual and with accessories other than those recommended by Antylia Scientific Ltd the protection provided might be impaired.

This equipment is designed to operate under the following conditions:

- For indoor use only
- Use in a well ventilated area
- Ambient temperature range +5°C to +40°C
- Altitude to 2000m
- Relative humidity not exceeding 80%
- Mains supply fluctuation not exceeding 10%
- Over-voltage category II IE C60364-4-443
- Pollution degree 2
- Use with a minimum distance all round of 200mm from walls or other items

The unit should be carried using both hands. Never move or carry the unit when in use or connected to the mains electricity supply.

Do not use in a hazardous atmosphere or with hazardous materials.

In the case of mains interruption, the unit will reset on restoration of the electricity supply.

### ELECTRICAL INSTALLATION

**This instrument must be earthed**

Before connection please read and understand these instructions and ensure that the line supply corresponds to that shown on the rating plate.

Voltage variants for these instruments are 230V 50Hz, 230V 60Hz, 120V 50Hz and 120V 60 Hz. The power consumption of each unit is 5W.

This unit is supplied with a mains lead which requires a rewirable mains plug to be fitted. The appropriate mains plug should be connected to the instrument **before** connection to the mains supply.

**It is important that this operation should only be undertaken by a qualified electrician.**

Note: Refer to the equipment's rating plate to ensure that the plug and fusing are suitable for the voltage and wattage stated.

The wires in the mains cable are coloured as follows:

| Connections | 220V-240V    | 110V-120V |
|-------------|--------------|-----------|
| Live        | Brown        | Black     |
| Neutral     | Blue         | White     |
| Earth       | Green/yellow | Green     |

Should the mains lead need replacement a cable of 0.5mm<sup>2</sup> to BS6500 & CENELEC HD21.5 should be used.

**IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN.**

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## SPECIFICATION

|  |  |  |
|--|--|--|
| Maximum time reading                                     | GT-5   | 9999 minutes in 1 minute increments    |
|  | GT-6   | 999.9 minutes in 0.1 minute increments |
| Amplitude of plunger oscillation                         | 13mm   |  |
| Frequency of plunger oscillation                         | GT-5   | 1 cycle/minute                         |
|  | GT-6   | 10 cycle/minute                        |
| Accuracy of gelation times                               | ±2%  |  |
| Overall dimensions, housing                              | H 119mm W 69mm D 94mm  |  |
| Overall length with standard plunger at bottom of stroke | 269mm  |  |
| Method of Support  | Tubular support 132mm long and 13mm diameter extends horizontally from rear of housing for attachment to retort stand. |  |
| Nett Weight  | 4.4kg  |  |
| Shipping Weight  | 7.7kg  |  |


## OPERATION

1. Mount the instrument clamping post onto a suitable retort stand or fixture at the desired height for operation. Ensure that the instrument is clamped firmly.
2. Screw the appropriate paddle to the end of the paddle connector on the sliding contact.
3. Connect the instrument to the electricity supply and switch on.
4. The instrument will perform an automatic self test at switch on. All segments will illuminate momentarily before displaying a single 'zero' on the right hand side., a buzzer will sound and the LED should illuminate green. The gelation timer is now ready to use.

## SAMPLE PREPARATION

The test requires about 100ml of sample in a container which has an internal diameter of at least 38mm. The top of the paddle stroke should be at least 32mm below the sample surface and the bottom of stroke not closer than 13mm to the bottom of the container. Standard sample cups for use with the instrument are available (catalogue number F7846).

## STARTING GELATION

Having prepared the sample in the manner described, gelation time measurement can be started in one of the following ways:

**If you want to set up your sample under the paddle, start reaction and gelation timer simultaneously:**

1. Press the 'start' button once. This starts the timer incrementing, the LED will flash 'Green' and the paddle will start to reciprocate (up and down). When the sample has gelled to the extent that **paddle movement is halted** (circuit is closed), the timer is stopped along with the paddle, the buzzer sounds and the display will alternate between the word 'GEL' and the indicated gelation time. This is the normal gelation mode.

**Alternatively, if you want to start the timer as soon as the sample reaction is started:**

2. Press the 'Hold' button once as soon as the sample is mixed, the timer starts incrementing (as gelation has started), the LED will illuminate 'Red', the paddle will reciprocate (up and down) and the display will alternate between the word 'HOLD' and accumulated time.

If the paddle switch is closed in the 'HOLD' state the timer is **not** stopped and the paddle continues to reciprocate. This allows you time to place the sample under the paddle and get ready for gelation measurement.

To commence gelation measurement once everything is set up, just press 'Start' button once and then follow instructions shown above in 1.

### **ACCIDENTAL TRIPPING DURING GELATION MEASUREMENT**

If the paddle is accidentally tripped during gelation, the timer can be restarted by pressing 'Restart' button once. Note: that the timer is restarted from when it was stopped, therefore the elapsed time between stopping and restarting the timer is not recorded.

### **STOPPING DURING GELATION MEASUREMENT**

The instrument can be stopped during gelation by pressing the 'Stop' button once. The paddle stops reciprocating, the LED is extinguished and the word 'STOP' is indicated on the display. To restart the gelation process just press the 'Restart' button once. Note: While in stop mode, the timer will continue to run, and when restarted will indicate actual gelation time. However, if the 'Stop' button is pressed a second time, having just stopped it, the timer will reset and the instrument will enter the 'Hold' state, (see 2. above).

## MAINTENANCE, SERVICING & REPAIR

**WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.**

This range of equipment does not require routine servicing. The only maintenance required is to clean external surfaces with a damp cloth and mild detergent. Spillages or splashes should be cleaned up straightaway after isolating the unit from the mains electricity supply.

### REPAIRS

Any repairs or replacement parts **must** be undertaken by suitably qualified personnel.

**Only spare parts supplied or specified by Antylia Scientific Ltd or its agent should be used. Fitting of non-approved parts may affect the performance of the safety features designed into the instrument.**

For a comprehensive list of parts required by service engineers conducting internal repairs, or to receive a service manual, please contact the Service Department of Antylia Scientific Ltd quoting both the model and serial number.

Antylia Scientific Ltd.

Service Department,

Beacon Road,

Stone,

Staffordshire,

ST15 0SA,

Email: [cpsservice@antylia.com](mailto:cpsservice@antylia.com).

Phone: +44 (0)1785 810475

If you require further technical or application assistance please contact:

Email: [cpstechsupport@antylia.com](mailto:cpstechsupport@antylia.com)

Phone: +44 (0)1785 810433

We are continually striving to improve our products and software. If you have any comments and suggestions on how we can do things better please send them to us at: [cpstechsupport@coleparmer.com](mailto:cpstechsupport@coleparmer.com)

### CALIBRATION

Your Gelation Timer has been carefully factory calibrated and is provided with a calibration certificate. If you require your unit to be re-calibrated at any time, please contact our Service Department of Antylia Scientific Ltd.

## ACCESSORIES

A range of accessories detailed below is available for use with the Gelation Timer.

| <b>Product No</b> | <b>Description</b>              |
|-------------------|---------------------------------|
| F0985             | 22mm Stainless steel plunger    |
| F0979             | 19mm Stainless steel plunger    |
| F0982             | 16mm Stainless steel plunger    |
| F1794             | Pack of 100 disposable plungers |
| F1795             | Disposable plunger adaptor      |
| F7846             | Pack of 240 sample cups         |

### DISPOSABLE PLUNGERS

In some applications the gel hardens so quickly that it is impossible to remove the plunger from the gel. In such work the use of Cole-Parmer disposable plungers is recommended. They are accurately weighed in accordance with Cole-Parmer specifications and made from glass. The plunger easily slides into the adaptor (F1795) screwed into the instrument in place of the stainless steel plunger.

#### Specification

|                 |           |        |
|-----------------|-----------|--------|
| Material        | Glass     |        |
| Disc dimensions | Diameter  | 21.8mm |
|                 | Thickness | 2.3mm  |
| Rod dimensions  | Diameter  | 4.5mm  |
|                 | Length    | 102mm  |

### STAINLESS STEEL PLUNGERS

For those applications where a small sample container is required, two alternative stainless steel plungers are available which are smaller than the standard plunger supplied with each instrument.

#### Specification

|                      | <b>Standard<br/>plunger<br/>F0985</b> | <b>Optional<br/>plunger<br/>F0982</b> | <b>Optional<br/>plunger<br/>F0979</b> |
|----------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| Disc dimensions (mm) |                                       |                                       |                                       |
| Diameter             | 22.4                                  | 16.0                                  | 19.1                                  |
| Thickness            | 2.3                                   | 2.5                                   | 2.3                                   |
| Rod dimensions (mm)  |                                       |                                       |                                       |
| Diameter             | 3.3                                   | 4.1                                   | 4.1                                   |
| Length               | 105.0                                 | 92.7                                  | 80.0                                  |
| Net weight (g)       | 11.4                                  | 11.4                                  | 11.4                                  |

## NOTES

This product meets the applicable CE Directives and UKCA Legislation for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in its vicinity will meet these standards and so we cannot guarantee

that interference will not occur in practise. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advise before use, contact the manufacturer.

Declaration of Conformity is available to view online at [www.coleparmer.com](http://www.coleparmer.com)

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T: +1 847 549 7600

**Cole-Parmer®**  
 essentials

Antylia Scientific Ltd  
 Beacon Road,  
 Stone,  
 Staffordshire,  
 ST15 0SA  
 United Kingdom