

SW-200 Series Slide Warmer

- Safely dry glass microscope slides after sample mounting in preparation for staining
- Holds up to 50 slides
- Slides can be arranged in different orientations
- Dark surface provides good visibility of samples
- 3-year warranty



SW-200 Series Slide Warmer

The SW-200 warming/drying bench aids in the preparation of microscope slides at the specimen mounting stage. It complements both the WD-200 wax dispenser and the FB-200 flotation bath to provide a complete set of paraffin wax sampling equipment for the modern laboratory.

Accepting up to 50 slides, the unit has the facility for drying slides in different orientations—across the drying support bars, angled from the bars, or flat without the bars. A silicone rubber mat element heater provides heating up to 100 °C (212 °F) and is controlled by a built-in energy regulator. The front panel has a temperature adjustment dial and two indicator lights for power and heater.



Specifications

Model no.	SW-200
Capacity	50 slides, 76 x 25 mm (3 x 1 in)
Temperature range	Ambient to 100 °C (212 °F)
Dimensions (W x H x D)	39.0 x 9.5 x 18.0 cm (15.4 x 3.7 x 7.1 in)
Warranty	3 years

Ordering Information

Description	Item no.	Series no.	Model no.	Legacy SKU
Microscope slide warmer, 230 VAC, 50/60 Hz	48573-06	SW-200	SW-200	MH6616
Microscope slide warmer, 115 VAC, 50/60 Hz	48573-01	SW-200	SW-200-115	MH6616X1

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1 800 363 5900
China: +86 21 5109 9909
France: +33 (0)1 48 63 78 00
Germany: +49 (0) 9377 92030

India: +1 800 266 1244
Italy: +39 (0)2 84349215
UK: +44 (0) 1480 272279
All other countries: +1 847 549 7600

Cole-Parmer®